


Prefab Orthoses

A cost effective option for at risk feet

This line of prefabricated insoles offers unique treatment solutions for a variety of control and protection needs often present in the Diabetic and Rheumatoid Arthritic populations. These prefab orthoses are available in three different designs; SOFT, GENTLE and MODERATE and can be an economical solution to provide interim protection and support while custom orthoses are being fabricated.

Paris

Orthotics



1 **Soft**



An off the shelf insole for those that require maximum cushioning and protection

- Smooth Poron top cover
- 15 durometer polyurethane insole
- Available in shoe sizes from Women's 5 to Men's 15

APPLICATIONS

- At risk feet that require exceptional cushioning and outstanding pressure relief
- A first response treatment prior to dispensing custom orthoses



2 **Gentle**



An off the shelf insole for those that require mild support, cushioning and protection

- P-Cell top cover
- 40 durometer EVA insole
- Available in shoe sizes from Women's 5 to Men's 15

APPLICATIONS

- To use when mild control and shear force diffusion is needed
- A first response treatment prior to dispensing custom orthoses



3 **Moderate**



An off the shelf insole for those that require support plus cushioning

- Woven nylon fabric and EVA top cover
- 55 durometer EVA insole
- Complementary Adhesive arch pads included
- Available in shoe sizes from Women's 5 to Men's 15

APPLICATIONS

- For individuals who need greater control, cushioning and insole durability
- Increased medial arch support when arch pads applied
- A first response treatment prior to dispensing custom orthoses

DISCLAIMER: Color of materials may change without notification





3630 East 1st Avenue
Vancouver, BC Canada
V5M 1C3

P 604.301.2150
T 800.848.0838
F 855.381.2157

parisorthotics.com

Prefab Orthoses



| SIZING | | | SOFT | | GENTLE | | MODERATE | |
|--|---------|---------|----------|--|----------|--|----------|--|
| Guide | Men's | Women's | Quantity | Price Per Pair | Quantity | Price Per Pair | Quantity | Price Per Pair |
| XS | 4 – 5 | 5 – 7 | | Please Note This order form is for Paris Orthotics Lab service customers ONLY . | | Retail Consumers: Please consult your clinician regarding where to purchase our prefab orthoses collection | | Lab Service Customer: Contact customer service or local account manager for pricing 1.800.848.0838 |
| SM | 6 – 7 | 8 – 9 | | | | | | |
| M | 8 – 9 | 10 – 11 | | | | | | |
| LG | 10 – 11 | 12 – 13 | | | | | | |
| XL | 12 – 13 | 14 + | | | | | | |
| XXL | 14 – 15 | | | | | | | |
| SUBTOTALS | | | # _____ | \$ _____ | # _____ | \$ _____ | # _____ | \$ _____ |
| TOTAL (Before Taxes & Shipping) | | | | | | | | \$ _____ |

VOLUME DISCOUNT – Orders over \$200 are eligible for a 10% discount

SHIP VIA

Include this order with my next orthotic shipment at no extra charge

Ship ASAP via UPS - bulk shipping charge will be added to your order

ORDER SUBMISSION OPTIONS

① **FAX** 1.855.381.2157

② **EMAIL** customerservice@parisorthotics.com

Request Confirmation:
 (please provide necessary contact information)

FAXED CONFIRMATION

EMAILED CONFIRMATION

SHIP TO

CLINICIAN _____

ACCOUNT No. _____

ACC. NAME _____

ADDRESS _____

CITY _____

PROVINCE _____ **POSTAL CODE** _____

PHONE _____

FAX _____

EMAIL _____