

# Prefab Orthoses



## A cost effective option for at risk feet

This line of prefabricated insoles offers unique treatment solutions for a variety of control and protection needs often present in the Diabetic and Rheumatoid Arthritic populations. These prefab orthoses are available in three different designs; SOFT, GENTLE and MODERATE and can be an economical solution to provide interim protection and support while custom orthoses are being fabricated.

Paris

Orthotics





- 1 **Soft**   | 1  
*An off the shelf insole for those that require maximum cushioning and protection*

- Smooth Poron top cover
- 15 durometer polyurethane insole
- Available in shoe sizes from Women's 5 to Men's 15

#### APPLICATIONS

- At risk feet that require exceptional cushioning and outstanding pressure relief
- A first response treatment prior to dispensing custom orthoses





- 2 **Gentle**   | 2  
*An off the shelf insole for those that require mild support, cushioning and protection*

- P-Cell top cover
- 40 durometer EVA insole
- Available in shoe sizes from Women's 5 to Men's 15

#### APPLICATIONS

- To use when mild control and shear force diffusion is needed
- A first response treatment prior to dispensing custom orthoses



- 3 **Moderate**   | 2 3  
*An off the shelf insole for those that require support plus cushioning*

- Woven nylon fabric and EVA top cover
- 55 durometer EVA insole
- Complementary Adhesive arch pads included
- Available in shoe sizes from Women's 5 to Men's 15

#### APPLICATIONS

- For individuals who need greater control, cushioning and insole durability
- Increased medial arch support when arch pads applied
- A first response treatment prior to dispensing custom orthoses

DISCLAIMER: Color of materials may change without notification





3630 East 1st Avenue  
Vancouver, BC Canada  
V5M 1C3

P 604.301.2150  
T 800.848.0838  
F 855.381.2157

parisorthotics.com

## Prefab Orthoses



SIZING			SOFT		GENTLE		MODERATE	
Guide	Men's	Women's	Quantity	Price Per Pair	Quantity	Price Per Pair	Quantity	Price Per Pair
<b>XS</b>	4 – 5	5 – 7		<b>Please Note</b> This order form is for Paris Orthotics Lab service customers <b>ONLY</b> .		<b>Retail Consumers:</b> Please consult your clinician regarding where to purchase our prefab orthoses collection		<b>Lab Service Customer:</b> Contact customer service or local account manager for pricing 1.800.848.0838
<b>SM</b>	6 – 7	8 – 9						
<b>M</b>	8 – 9	10 – 11						
<b>LG</b>	10 – 11	12 – 13						
<b>XL</b>	12 – 13	14 +						
<b>XXL</b>	14 – 15							
SUBTOTALS			# _____	\$ _____	# _____	\$ _____	# _____	\$ _____
<b>TOTAL</b> (Before Taxes & Shipping)								\$ _____

**VOLUME DISCOUNT** – Orders over \$200 are eligible for a 10% discount

**SHIP VIA**

Include this order with my next orthotic shipment at no extra charge

Ship ASAP via UPS - bulk shipping charge will be added to your order

**ORDER SUBMISSION OPTIONS**

① **FAX** 1.855.381.2157

② **EMAIL** customerservice@parisorthotics.com

**Request Confirmation:**  
 (please provide necessary contact information)

**FAXED CONFIRMATION**

**EMAILED CONFIRMATION**

**SHIP TO**

**CLINICIAN** \_\_\_\_\_

**ACCOUNT No.** \_\_\_\_\_

**ACC. NAME** \_\_\_\_\_

**ADDRESS** \_\_\_\_\_

\_\_\_\_\_

**CITY** \_\_\_\_\_

**PROVINCE** \_\_\_\_\_ **POSTAL CODE** \_\_\_\_\_

**PHONE** \_\_\_\_\_

**FAX** \_\_\_\_\_

**EMAIL** \_\_\_\_\_