

Carbon Fibre Spring Plate

Internal Shoe Stiffener

Paris

A simple and effective solution for limiting motion across the forefoot and aiding in propulsion



FEATURES:

- Extremely thin, light and durable made from premium quality pre-preg carbon
- Pre-shaped left or right and available in 4 sizes
- Contoured design to provide ideal toe spring and heel elevation for most athletic, walking and orthopaedic footwear designs
- Excellent interface with custom foot orthotics
- A practical alternative to some common shoe modifications

CLINICAL INDICATIONS:

- Post-surgical care option when transitioning from post-op shoe to regular footwear after bunion, hammertoe and other forefoot surgeries.
- Management solution for turf toe, forefoot stress fractures, sesamoiditis and Lisfranc injuries.
- Protection and toe off assistance for pathologies affecting the midfoot and forefoot including rheumatoid arthritis, hallux rigidus or limitus, Freiberg's disease, diabetes and forefoot amputations.

Sizing: Plates are available left and right in 4 sizes

Orthotics



CARBON FIBRE SPRING PLATE

CONSIDERATIONS:

- Maximum body weight tolerance is 280 lbs.
- Important to use in shoes with adequate forefoot rocker profile.
- Plates are not returnable for credit however size exchange requests will be accommodated for **UNUSED AND UNMODIFIED PLATES ONLY.**

PRODUCT GUARANTEE

The Paris Carbon Fibre Spring Plate is guaranteed against breakage for 6 months.

- Minor grinding (forefoot shape and/or length) will not affect durability or nullify warranty.
- Heating and reshaping will cause delaminating and will nullify warranty.

SHIP TO

ACC. NAME _____

ACCOUNT # _____

CLINICIAN _____

ADDRESS _____

CITY _____

PROVINCE _____ POSTAL CODE _____

PHONE _____ FAX _____

EMAIL _____

Bill To Orders placed without an Account # are payable by credit card. Provide payment information below or call customer service to pay-by-phone. **1.800.848.0838**

CARD No. _____ - _____ - _____ - _____

CARD TYPE VISA MC AMEX EXPIRY ____ / ____

ORDER SUBMISSION OPTIONS

- ① FAX - 1.855.381.2157
- ② EMAIL - customerservice@parisorthotics.com

PLEASE NOTE:
 All received faxes will be sent a confirmation

SHIP VIA

Include this order with my next orthotic shipment at no extra charge

Ship ASAP via UPS - bulk shipping charge will be added to your order

ORDER DETAILS – SIZING & UNITS

SIZE	DEVICE LENGTH tolerance +/- 2MM	INSOLE MEASUREMENT	PRICE / UNIT	*BULK ORDER PRICE / UNIT	LEFT / RIGHT	QTY	ITEM TOTAL (\$)	
XS	239MM					L		
S	259MM					L		
M	279MM					L		
L	298MM					L		
XS	239MM					R		
S	259MM					R		
M	279MM					R		
L	298MM					R		

* Orders of 8 or more units are eligible for a \$5.00 discount per unit.
 Recommended **STARTER KIT**: one of each size for L & R.

SUB-TOTAL
 plus applicable taxes