

INTRODUCING THE



GAUNTLET ANKLE FOOT ORTHOSES

A FUNCTIONAL GAUNTLET FROM A NEUTRAL SUSPENSION CAST

FEATURES:

- Medial and Lateral Arch Suspender* Legitimate Varus / Valgus control of the hind foot
- Neutral Suspension Casts Accepted
- Non-weight bearing cast provides better heel and arch contour
- All casts are intrinsically balanced to correct forefoot varus / valgus deformities

*US. Patent No. 6,602,215

LEGITIMATE VARUS / VALGUS CONTROL OF THE HIND FOOT

CLINICAL INDICATIONS:

- Rigid, non-reducible Adult Acquired Flatfoot (Stage III & IV)
- Severe Degenerative Joint Disease or Deformity of Hind Foot
- Charcot Arthropathy

FEATURES:

- Varus or Valgus lifting strap
- Intrinsic balanced Ankle Foot Orthoses (AFO) shell
- Extrinsic rearfoot post

Use of the finest leather and extra padding assures patient comfort and support



Padded top-collar to reduce calf irritation

Custom molded AFO shell for superb comfort and support

Arch cut out with lifting strap provides adjustable control of hind foot

Specially designed Heel Post adds stability



Medial Arch Suspender (9")



Lateral Arch Suspender (7")

ORDERING OPTIONS:

- Tan – 7" height
- Tan – 9" height
- Chocolate – 7" height
- Chocolate – 9" height



Tan** (7")




Chocolate (9")

**The tan Gauntlet images showcased here are for marketing purposes only. They feature non-standard contrasting elements chosen specifically to emphasize the various features of the brace. For ordering purposes, please note that the lab standard tan gauntlet comes with corresponding tan cross strap and tan laces. **

PLEASE NOTE: The Richie Gauntlet is a SPECIAL ORDER item. Shipping turn around times average 10-14 BUSINESS DAYS.

THE **Richie GAUNTLET™**


THE RICHIE GAUNTLET AFO PRESCRIPTION FORM

 <p>Paris Orthotics Ltd. Head Office / Lab Division 3630 E. 1st Ave Vancouver BC Canada V5M 1C3 P 604.301.2150 F 604.301.2157</p>	Clinician:
	Account:
	Address:
	City: Province: Postal:
	PO: (internal use ONLY)
	Pt Name: <input type="checkbox"/> Male <input type="checkbox"/> Female Age:
	Height: Weight: Shoe Size:
	Cast enclosed for <input type="checkbox"/> Left <input type="checkbox"/> Right <input type="checkbox"/> B/L
	<ul style="list-style-type: none"> • please landmark MEDIAL & LATERAL MALLEOLI on the negative cast • cast must have a standard height <u>GREATER than 9"</u>

CLINICAL INFORMATION

DIAGNOSIS:	Location(s) Accommodation (describe & landmark on cast):
-------------------	---

PRESCRIPTION

<p>Height:</p> <p><input type="checkbox"/> 7" – most versatile height</p> <p><input type="checkbox"/> 9" – For maximal rigidity and control</p> <p>ARCH SUSPENDER: (SELECTION REQUIRED)</p> <p><input type="checkbox"/> Medial (varus force on hindfoot)</p> <p><input type="checkbox"/> Lateral (valgus force on hindfoot)</p> <p><input type="checkbox"/> None</p>	<p>Color:</p> <p><input type="checkbox"/> Tan</p> <p><input type="checkbox"/> Chocolate</p>	 <p>7" (left) 9" (right)</p>
--	--	---

CASTING INSTRUCTIONS USING THE STS MID LEG SOCK



Cutting strip, tubing & bag in place



Gather STS Sock-slide onto foot/ leg



Position STJ neutral: Ankle at 90°



OR: semi weight bearing on foam



Cut along cutting strip-full length



Mark medial & lateral malleolus