

The Strassburg Sock

User-Friendly Night Splint

Paris

SIZE OPTIONS

- Adult Regular (up to 16" calf diameter)
- Adult Large (>16" calf diameter)
- A Pediatric version is also available for use in treating ankle equinus

The Strassburg sock is a very effective, surprisingly comfortable and cost effective compliment to traditional plantar fasciitis treatment regimes.

Conventional hard plastic night splints maintain tension on the fascia by holding the ankle at or near a neutral position. The Strassburg sock provides a more complete and specific tension on the plantar fascia by passively dorsiflexing the toes (and engaging the windlass mechanism) as well as by holding the ankle in neutral position. Used in combination with custom foot orthotic, anti-inflammatory and other conservative therapies, the Strassburg sock night splint significantly improves clinical outcomes.

THE STRASSBURG SOCK ADVANTAGE:

- A breathable and washable tubular knit "sock" design that is comfortable enough to sleep in
- A unique hook & loop adjustable strap design that dorsiflexes the toes as well as holds the foot in a neutral position
- No restriction of the ankle joint



A retrospective study published by Barry and colleagues¹, compared the Strassburg sock to a gastrocnemius-soleus stretching regime reported the following results:

	STRASSBURG SOCK GROUP	STRETCHING REGIME GROUP
# OF SUBJECTS	89 patients	70 patients
MEAN RECOVERY TIME	18.5 days	58.6 days
FOLLOW UP VISITS	1.78	3.07
ADDITIONAL INTERVENTIONS	1.83	2.14

Of the 89 patients in the Strassburg sock group, all tolerated wearing the intervention and 97.8% of them recovered within 8 weeks.

¹ Barry LD et al. A Retrospective Study of Standing Gastrocnemius-soleus Stretching versus Night Splinting in the Treatment of Plantar Fasciitis. *The Journal of Foot & Ankle Surgery*. 41 (4):221-227, 2002.

Orthotics

THE STRASSBURG SOCK



***The Pediatric Sock is a SPECIAL ORDER item**
 Special order shipping turnaround times average 10-14 working days. All other sizes are readily available.

SHIP TO

CLINICIAN _____

ACCOUNT No. _____

ACC. NAME _____

ADDRESS _____

CITY _____

PROVINCE _____ POSTAL CODE _____

PHONE _____

FAX _____

EMAIL _____

ORDER SUBMISSION OPTIONS

- ① FAX - 604.301.2157*
 - ② EMAIL - customerservice@parisorthotics.com
- PLEASE NOTE:**
 *All received faxes will be sent a confirmation

SHIP VIA

Include this order with my next orthotic shipment at no extra charge

Ship ASAP via UPS - bulk shipping charge will be added to your order

ORDER DETAILS – SIZING & UNITS

QTY	SIZE	\$ / UNIT	ITEM TOTAL
	Regular (up to 16" calf diameter)	This order form is for Paris Orthotics Lab service customers ONLY.	
	Large (>16" calf diameter)		
	Pediatric (up to 12" calf diameter) SPECIAL ORDER*		

VOLUME DISCOUNT – Price Breaks

- 6 – 12 (assorted variety) **5% Discount**
- 13 – 18 (assorted variety) **10% Discount**
- 19+ (assorted variety) **15% Discount**

DISCOUNT

SUB-TOTAL

plus applicable taxes